

## Flat Pouches

The wide range of standard sizes allows optimum choice of a correct-sized pack for each item and for double packing purposes. Pouches are printed with process indicators for Steam, EO and Formaldehyde sterilization.



Article	Width	Length	Packing Unit
S1	50	250	2400
S10	205	400	600
S10A	205	440	600
S11	250	500	600
S12	100	570	1800
S13	150	400	1200
S14	300	570	600
S15	150	200	1200
S16	250	380	1200
S17	50	200	3600
S18	75	300	4200
S19	100	350	1200
S2	75	200	2400
S20	150	300	1800
S2070	205	700	600
S21	300	500	600
S22	75	520	2400
S2244	220	440	600
S23	75	230	1800

Article	Width	Length	Packing Unit
S24	75	150	2400
S25	100	150	2400
S26	150	350	1800
S27	120	400	1200
S28	150	520	1200
S29	270	350	600
S3	100	200	1800
S30	270	440	600
S31	420	500	500
S32	420	600	500
S33	180	320	1000
S34	160	600	1200
S35	205	600	1200
S36	250	600	500
S37	250	700	600
S38	320	500	600
S4	75	270	3600
S5	100	270	3000
S6	150	270	1800
S7	205	270	600
S7A	205	340	1000
S8	100	300	1200
S9	100	400	1200