

Tray Liners



Tray Liners for use in surgical instrument trays for absorbing moisture and for mechanical protection

- Absorbs excessive condensation.
- Wet strengthened special paper.
- Virtually lint-free.

TL - standard Trayliners - into a tray under the instruments

TLWE - Trayliners with edge - outside of the instrument tray on
the top of the inner wrap

TLNW - Nonwoven Trayliners - for combinations of nonwovens
and SMS-type of wraps - extreme absorption capacity

Article	g/m ²	Width	Length	Packing Unit
TL25	70	250	300	1000
TL3040	70	300	400	
TL40	70	400	500	500
TL50	70	300	500	500
TLNW1525	45	150	250	3200
TLNW3025	45	300	250	1200
TLNW3040	45	300	400	1200
TLNW3050	45	300	500	600

Article	g/m2	Width	Length	Packing Unit
TLWE2525	107	250	250	500
TLWE4825	107	480	250	500
TLWE5030	107	500	300	500