

Sterilization Tapes



- ITSG and JTBSG with STEAM green-to-brown indicator. In natural and blue base material. ITSG also available as latex free version called ITSGL.
- ITE with EO gas purple-to-green indicator.
- ITH with HOT AIR green-to-brown indicator.
- ITP tape strips with VH2O2 blue-to-pink indicator.
- ITP25 tape reel with VH2O2 blue-to-pink indicator.
- Neutral closing tapes with no imprint. In natural and blue base material.

All Indicators are process indicators in accordance with ISO 11140-1 Type 1.

Article	Width	Length	Gusset	Packing Unit
ITE19	19			48
ITH19	19			48
ITN19	19			48
ITN25	25			36
ITN50	50			24
ITP20	20	150		200
ITP25	25			24
ITS12G	12			72
ITS12GL	12			72
ITS19G	19			48
ITS19GL	19			48
ITS25G	25			36

Article	Width	Length	Gusset	Packing Unit
ITS25GL	25			36
JTB19N	19			48
JTB19SG	19			48
JTB25N	25			36
JTB25SG	25			36
TDS25				1