

Sterilization Tapes



- ITSG and JTBSG with STEAM green-to-brown indicator. In natural and blue base material. ITSG also available as latex free version called ITSGL.
- ITE with EO gas purple-to-green indicator.
- ITH with HOT AIR green-to-brown indicator.
- ITP tape strips with VH2O2 blue-to-pink indicator.
- ITP25 tape reel with VH2O2 blue-to-pink indicator.
- Neutral closing tapes with no imprint. In natural and blue base material.

All Indicators are process indicators in accordance with ISO 11140-1 Type 1.

Article	Width	Length	rm/roll	Packing Unit
ITE19	19		50	48
ITH19	19		50	48
ITN19	19		50	48
ITN25	25		50	36
ITN50	50		50	24
ITP20	20	150		200

Article	Width	Length	rm/roll	Packing Unit
ITP25	25		25	24
ITS12G	12		50	72
ITS12GL	12		50	72
ITS19G	19		50	48
ITS19GL	19		50	48
ITS25G	25		50	36
ITS25GL	25		50	36
JTB19N	19		50	48
JTB19SG	19		50	48
JTB25N	25		50	36
JTB25SG	25		50	36
TDS25				1