

STERIKING®

WIPAK

Health

Steriking® LT-Blueline

Pouches and Rolls for Low Temperature sterilization



High protection

- Strongest medical style Tyvek® 1073B (75g/m²) and 9 layer film for extra protection

Outstanding resistance for microbial penetration

- Superior bacterial barrier and protection against contamination with moisture resistance

Efficient sterilization

- Allows quick passage and escape of sterilizing agent into/from the package


Aseptic presentation

- Clean fiber free opening

Steriking® LT-Blueline

Developed for packing items for low temperature sterilization

Comparison of technical properties




Suitable for gas (EO, FO), hydrogen peroxide and irradiation sterilization

- Validation studies for compatibility with intended sterilization processes
- Shelf life studies
- Double pouching studies

Ref.	Description
LTS	Heat-Sealable Flat pouches
LTSS	Self-Sealable Flat pouches
LTR	Flat rolls in 100m length



Reference Color Chart for LT-Blueline



We care that you pack safely!



Involve. Inspire. Impress.

www.steriking.com
e-mail: steriking@wipak.com