

Sterilization Tapes



- ITSG and JTBSG with STEAM green-to-brown indicator.
In natural and blue base material. ITSG also available as latex free version called ITSGL.
- ITH with HOT AIR green-to-brown indicator.
- ITP tape strips with VH2O2 blue-to-pink indicator.
- Neutral closing tapes with no imprint. In natural and blue base material.

All Indicators are process indicators in accordance with ISO 11140-1 Type 1.

Article	Width	Length	Gusset	Packing Unit
ITH19	19			48
ITN19	19			48
ITN25	25			36
ITN50	50			24
ITP20	20	150		200
ITS19G	19			48

Article	Width	Length	Gusset	Packing Unit
ITS19GL	19			48
ITS25G	25			36
ITS25GL	25			36
JTB19N	19			48
JTB19SG	19			48
JTB25N	25			36
JTB25SG	25			36
TDS25				1