

# Sterilization Tests



## For routine control of autoclave performance:

- Indicator sheets, Daily Test Sheets, for traditional Bowie & Dick type test with a cloth package in accordance with ISO 11140-1 (type 2) and ISO 11140-3 for pre-vacuum steam sterilizers to demonstrate compliance with EN 285. Colour changing from white to black.
- Ready made, compact, single-use Bowie & Dick type test packs: DC20 and DCG2.
- Helix Challenge Test as a product challenge device (PCD) for regular control of hollow loads in autoclaves. Helix consists of the device and 100 indicators strips.  
Compliance with ISO 11140-1 (type 6) and EN 867-5.

Article	Width	Length	Gusset	Packing Unit
DC200	135	105	20	20
DCG2	100	150	10	20
DT50	210			50
HX100	115	297	23	25
HX400	115	175	23	25